**ergodyne®**

INTERGALACTIC HEADQUARTERS  
1021 BANDANA BOULEVARD EAST // SUITE 220 // SAINT PAUL, MN 55108  
TOLL-FREE 800 225 8238 // INT'L +1 651 642 9889 // ERGODYNE.COM

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
## REUSABLE HOT/COLD PACK

### COLD APPLICATION

To help reduce swelling, inflammation and muscle spasms after a sprain, strain or minor injury.

- Place pack flat in freezer for at least 2 hours.
- Place pack in a towel or into the support product before application to protect skin. Never use the pack directly on skin.
- Leave cold pack on for 15-20 minutes. Return pack to freezer for 30 minutes between uses if necessary. Repeat every 1 – 1 ½ hours as needed. Consult a physician about extended use.

Never apply cold pack for more than 30 minutes at a time. Direct or extended application of cold pack can cause frostbite. Check skin for redness or discoloration during use, if skin appears red or discolored stop treatment at once.

### HEAT APPLICATION

Targets muscle pain and tension, stiffness and chronic pain.

Note: Heated pack should be warm but not hot and never painful to touch. If the pack is too hot let it cool at room temperature for several minutes and check it again before application.

Do not use a damaged, ruptured, or leaking pack. Allow a damaged pack to cool down in an inaccessible location and then dispose of it in the garbage. Do not wash gel down drains.

#### HOT WATER METHOD

Boil approximately 4 -5 liters (4-5 quarts) of water.

- Remove water from heat.
- Immerse room temperature pack in water for approx. 10 minutes.
- Remove pack from water using tongs. Never use sharp utensils.
- Dry the pack and check for desired temperature. Let cool if necessary.
- Place pack in a towel or into the support product before application to protect skin. Never use the pack directly on skin.
- Apply in 20-minute intervals. Consult a physician regarding extended use.

#### MICROWAVE METHOD

Note: Microwaves vary in power and may heat unevenly. If the pack is heated unevenly, knead the gel or gently shake the pack to better distribute heat inside the pack.

If the pack expands like a balloon stop the microwave immediately. The pack has become overheated. Do not use an overheated pack, let the pack cool and discard it.

- Pack must be at room temperature before use. Flatten to ensure gel is evenly distributed.
- Heat on High for 50 secs and let stand for 50 secs. If the pack is not warm continue heating in 10 second intervals until a safe and desired temperature is reached.
- Place pack in a towel or into the support product before application to protect skin. Never use the pack directly on skin.
- Apply in 20-minute intervals. Consult a physician regarding extended use.

#### CAUTION:

DO NOT APPLY DIRECTLY TO SKIN.

Without protection skin may become raw or damaged from the cold, or severely burned from heat. Check the temperature of the pack before using. If it is too hot or cold let it sit at room temperature for several minutes before checking again. ALWAYS USE A TOWEL TO COVER PACK OR INSERT IT INTO THE PROVIDED SUPPORT PRODUCT BEFORE USE.

Do not lean, sit, or lie on the pack as it may break and leak. Avoid contact with gel. Ruptured packs and gel should be placed in the garbage, do not wash gel down drains. Clothing and support products should be washed thoroughly with soap and water to remove gel residue. Wash hands thoroughly with soap and water after contact with gel.

Heat and cold packs should not be used by individuals who are sensitive to temperature extremes, such as those with circulatory problems, diabetes, nerve damage, paralysis, or sensitive skin.

#### USE ONLY AS DIRECTED.

KEEP OUT OF REACH OF CHILDREN.

### INSERTING HOT/COLD PACK INTO SUPPORT

- Follow the instructions on the Hot/Cold Pack to activate it for either warming or cooling (Fig. 1)
- Slide the Hot/Cold Pack into the mesh pocket on the inside of the Back Support so that the slotted tabs are at the top. (Fig. 2)
- Use the buttons to secure the Hot/Cold Pack into the mesh pocket via the slotted tabs. (Fig. 3)
- Your Back Support is now ready to provide warming or cooling relief. (Fig. 4)

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©2023 ERGODYNE. ALL RIGHTS RESERVED. ALL WRONGS REVERSED.